



Garment Measuring Guide

HUNTLEIGH


A MEMBER OF THE ARJO FAMILY

Garment Measuring Guide

The following guide is intended as an aid to help the user select the correct size garment.

When determining the correct size garment:

- Remove any bandaging, stockings or bulky dressings
- Measure the points as shown on the diagram below
- DO NOT pull the measuring tape tight
- Make sure the garment extends beyond the area under treatment.



*Measurements should be made at the widest point on the calf and thigh.

Lymph Assist (Hydroven 12) Arm

Twelve Chamber Arm Garments



Order Code	Type	Length (L)	Shoulder Circumference (S)	Hand Circumference (H)
1533523	Full Arm	68cm	62cm	44cm
1533524	Full Arm	78cm	62cm	44cm

Lymph Assist (Hydroven 12) Arm Insert

Twelve Chamber Arm Inserts



Order Code	Type	Length (L)	Additional Width at Shoulder	Additional Width at Hand
1533527	Full Arm	68cm	17cm	13cm
1533528	Full Arm	78cm	17cm	13cm

Lymph Assist (Hydroven 12) Leg

Twelve Chamber Leg Garments



Order Code	Type	Length (L)	Thigh Circumference (T)	Calf Circumference (C)
1533519	Standard Leg	76cm	71cm	59cm
1533520	Standard Leg	84cm	71cm	59cm
1533521	Wide Leg	76cm	79cm	62cm
1533522	Wide Leg	84cm	79cm	62cm

Lymph Assist (Hydroven 12) Leg Insert

Twelve Chamber Leg Inserts



Order Code	Type	Length (L)	Additional Width at Thigh	Additional Width at Foot
1533525	Leg Insert	76cm	19cm	14cm
1533526	Leg Insert	84cm	19cm	14cm



2019