



Keeping Health Simple™

AllBrand™ ECG Electrodes Belt


12-lead | Connect to All ECG Brands | Fast & Easy to Place
One Size Fits All - Elastic | For Men & Women | Reusable
No Shaving, No Pain, No Disposables



AllBrand™ ECG Electrodes Belt


Compatible with **All 12-lead ECG Brands**

AllBrand™ with Handles



Ref 11005/h

AllBrand™ with Straps



Ref 11005/s

- ✓ One size fits all - Elastic
- ✓ Intended for 12-lead Rest ECG
- ✓ Compatible with all 12-lead ECG brands
- ✓ Connect to lead wires using Snap or Clip Connectors
- ✓ Built-in 6 chest Ag/AgCL electrodes
- ✓ With Handles / Straps
- ✓ No maintenance
- ✓ Electrodes and Plastic Parts may be replaced on-site
- ✓ With Hygienic protection

- ✓ Provided with:

- 4 limb clamp Ag/AgCl electrodes
- Banana to Clip adaptors (for 3&4 mm leadwires)
- Soft Carry Case



- ✓ 12 month warranty

Warranty provided for silicone parts.



* Note: All AllBrand™ ECG Belts are supplied without ECG patient cable.